

Anti-Ferritin Heavy Chain Rabbit antibody



Catalog Number: 169064

Product name

Polyclonal

Anti-Ferritin Heavy Chain Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

ELISA

Antibody description

Ferritin Heavy Chain Rabbit polyclonal antibody

Validations

Preparation

Antigen: Recombinant protein of human Ferritin.


Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality



Standard Curve for Ferritin: Capture Antibody Mouse mAb (201067, Ferritin Light Chain Mouse mAb) to Ferritin at 2u03bcg/ml and Detector Antibody Rabbit pAb(3010667, Ferritin rabbit pAb) to Ferritin at 0.5u03bcg/ml.