

# CD2 antibody

Catalog Number: 109064

## Product name

CD2 antibody

## Specificity

Human; other species not tested.

## Antibody description

CD2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IHC detected in human tonsillitis tissue. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of CD2 recombinant protein (Accession Number: NM\_001767). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109064(CD2 Antibody) at dilution of 1:3000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109064(CD2 Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieved with Tris-EDTA buffer, pH9.0



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109064(CD2 Antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieved

# CD2 antibody



Catalog Number: 109064

---

with Tris-EDTA buffer, pH9.0