

HRH2 antibody



Catalog Number: 111464

Product name

HRH2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HRH2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, human skeletal muscle tissue, mouse stomach tissue. Observed molecular weight by Western-blot: 59kd,69kd

Preparation

This antibody was obtained by immunization of HRH2 recombinant protein (Accession Number: NM_001131055). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:111464(HRH2 antibody) at dilution of 1:500