

# RAB9A-Specific antibody



Catalog Number: 114464

## Product name

RAB9A-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RAB9A-Specific Rabbit Polyclonal antibody.  
Positive WB detected in HEK-293 cells, HeLa cells.  
Positive IHC detected in human gliomas tissue.  
Observed molecular weight by Western-blot: 23 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: XM\_047442644).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114464(RAB9A-Specific antibody) at dilution of



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:114464(RAB9A-Specific Antibody) at dilution of 1:50