

Anti-CLEC6A antibody



Catalog Number: 176464

Product name

Anti-CLEC6A antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IHC-P

Antibody description

Rabbit polyclonal antibody to CLEC6A

Dilutions

WB: 1:500-1:2000

Preparation

This antigen of this antibody was synthesized peptide derived from the internal region of human dectin-2.

IHC-P: 1:100-1:300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations




Fig1:: Western blot analysis of SH-SY5Y 293T using Dectin-2 antibody. Antibody was diluted at 1:1000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000