

# Anti-Scavenging Receptor SR-BI Rabbit antibody



Catalog Number: 166864

## Product name

Rabbit IgG

Anti-Scavenging Receptor SR-BI Rabbit antibody

## Applications

## Specificity

WB, ICC/IF, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

Scavenging Receptor SR-BI Rabbit polyclonal antibody

ICC/IF: 1/20

## Preparation


IP: 1/20

Antigen: Recombinant protein of human Scavenging Receptor SR-BI

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Western blot detection of Scavenging Receptor SR-BI in 3T3, HeLa cell lysates using Scavenging Receptor SR-BI Rabbit pAb(1:1000 diluted). Predicted band size: 61 kDa. Observed band size: 80 kDa.