

# ETV6 Polyclonal Antibody



Catalog Number: 163164

## Product name

ETV6 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ETV6

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ETV6 (NP\_001978.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IP


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IP 1:50 - 1:100

## Validations



Western blot - ETV6 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ETV6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 90s.

Immunoprecipitation - ETV6 Polyclonal Antibody

Immunoprecipitation analysis of 200ug extracts of SW620 cells using 3ug ETV6 antibody . Western blot was performed from the immunoprecipitate using ETV6 antibody at a dilution of 1:500.