

# CLTCL1 antibody



Catalog Number: 107164

## Product name

CLTCL1 antibody

## Specificity

Human; other species not tested.

## Antibody description

CLTCL1 Mouse Monoclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells. Observed molecular weight by Western-blot: 180 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_007098). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with (CLTCL1 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:107164(CLTCCL1 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)