

FBXW4 antibody



Catalog Number: 110564

Product name

FBXW4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FBXW4 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, mouse liver tissue. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of FBXW4 recombinant protein (Accession Number: BC007380). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110564(FBXW4 Antibody) at dilution of 1:300