

# SPERT antibody



Catalog Number: 115564

## Product name

SPERT antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SPERT Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in mouse testis tissue, mouse spermatophore tissue, PC-3 cells. Observed molecular weight by Western-blot: 48-52 kDa

## Preparation

This antibody was obtained by immunization of SPERT recombinant protein (Accession Number: NM\_001286342). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115564(SPERT antibody) at dilution of 1:1000



Immunofluorescent analysis of Hela cells, using SPERT antibody Catalog No:115564 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).