

UHRF1BP1L antibody



Catalog Number: 116564

Product name

UHRF1BP1L antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

UHRF1BP1L Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 160-165 kDa

Preparation

This antibody was obtained by immunization of UHRF1BP1L recombinant protein (Accession Number: BC127017). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116564(UHRF1BP1L Antibody) at dilution of



1:1500

IP Result of anti-UHRF1BP1L (IP:Catalog No:116564, 5ug; Detection:Catalog No:116564 1:1000) with mouse brain tissue lysate 5000ug.