

# PLEKHO1 antibody



Catalog Number: 113964

## Product name

PLEKHO1 antibody

## Immunogen

[Human PLEKHO1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PLEKHO1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse testis tissue. Observed molecular weight by Western-blot: 46 kDa

## Preparation

This antibody was obtained by immunization of PLEKHO1 recombinant protein (Accession Number: NM\_016274). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113964(PLEKHO1 Antibody) at dilution of 1:300