

HB9 antibody



Catalog Number: 111264

Product name

Polyclonal

HB9 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

HB9 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue, mouse liver tissue. Observed molecular weight by Western-blot: 45-50 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005515). Purification method: Antigen affinity purified.

Formulation


PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111264(HB9 Antibody) at dilution of 1:600