

# LACTB antibody



Catalog Number: 112264

## Product name

LACTB antibody

WB: 1:1000-1:10000

IHC: 1:20-1:200

## Immunogen

[Human LACTB Recombinant protein \(His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

LACTB Rabbit Polyclonal antibody. Positive IHC detected in human liver. Positive WB detected in human heart tissue, HeLa cells. Observed molecular weight by Western-blot: 55kd

## Preparation

This antibody was obtained by immunization of LACTB recombinant protein (Accession Number: NM\_032857). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112264(LACTB antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human liver using Catalog No:112264(LACTB antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:112264(LACTB

# LACTB antibody



Catalog Number: 112264

---

antibody) at dilution of 1:100 (under 40x lens)