

# TPH2 antibody

Catalog Number: 116264

---

## Product name

TPH2 antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

TPH2 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 56 kDa

## Preparation

This antibody was obtained by immunization of TPH2 recombinant protein (Accession Number: BC114442). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116264(TPH2 Antibody) at dilution of 1:300