

Anti-Nectin 2 Rabbit antibody



Catalog Number: 167264

Product name

Anti-Nectin 2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Nectin 2 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human Nectin 2

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP, FC

Dilutions


WB: 1/500-1/2000

ICC/IF: 1/50-1/200

IP: 1/100

FC: 1/50

Validations



Western blot analysis of Nectin 2 expression in SW480 cell lysate.