

HP antibody



Catalog Number: 107264

Product name

HP antibody

Specificity

Human, Pig; other species not tested.

Antibody description

HP Mouse Monoclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of HP recombinant protein (Accession Number: BC058031). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107264(HP Antibody) at dilution of 1:1000