

Anti-Syntenin (3D9) Mouse antibody



Catalog Number: 168264

Product name

Anti-Syntenin (3D9) Mouse antibody

Specificity

Human

Antibody description

Syntenin (3D9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Syntenin protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Syntenin in HeLaNE and HeLa cell lysates using Syntenin mouse mAb (1:1000 diluted). Predicted band size: 32KDa. Observed band size: 32KDa.