

# Anti-Syntenin (3D9) Mouse antibody



Catalog Number: 168264

## Product name

Anti-Syntenin (3D9) Mouse antibody

## Ig Type

Mouse IgG2b

## Specificity

Human

## Applications

WB

## Antibody description

Syntenin (3D9) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant human Syntenin protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot detection of Syntenin in HeLaNE and HeLa cell lysates using Syntenin mouse mAb (1:1000 diluted). Predicted band size: 32KDa. Observed band size: 32KDa.