

Anti-PTPRB antibody



Catalog Number: 176664

Product name

Anti-PTPRB antibody

Ig Type

Mouse IgG2a

Specificity

Human, Mouse

Applications

ICC

Antibody description

Mouse monoclonal antibody to PTPRB

Dilutions

ICC: 1:1,000

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Validations




Fig1: ICC staining PTPRB in F9 cells (red). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS and counterstained with DAPI in order to highlight the nucleus (blue).