

# Anti-ABCB5 antibody



Catalog Number: 100364

## Product name

Anti-ABCB5 antibody

## Specificity

Human ABCB5

## Antibody description

Rabbit polyclonal to ABCB5

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the human ABCB5, and purified by antigen affinity chromatography.

## Formulation

0.2  $\mu$ m filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 10-20  $\mu$ g/ml

## Validations



Items	Lanes	A
Sample (Membrane Lysate)		K562
Sample Volume ( $\mu$ g/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		10-20 $\mu$ g/ml
Secondary Antibody		Dylight 800 labeled Antibody To Rabbit IgG (H+L) at 1:5000 dilution
Developed using Odyssey imaging system.		
Explanation		Predicted band size : 90 kDa Observed band size : 90 kDa

ABCB5 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot