

# HIST3H2A antibody



Catalog Number: 111364

## Product name

HIST3H2A antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

HIST3H2A Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 14-18kd

## Preparation

This antibody was obtained by immunization of HIST3H2A recombinant protein (Accession Number: BC001193). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111364(HIST3H2A antibody) at dilution of 1:100