

# Anti-NLGN3 antibody



Catalog Number: 142364

## Product name

Anti-NLGN3 antibody

## Immunogen

[Human Neuroligin 3 / NLGN3 \(His Tag\) recombinant protein](#)

## Specificity

Human Neuroligin 3 / NLGN3

## Antibody description

Rabbit Monoclonal to Human NLGN3

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Neuroligin 3 / NLGN3 .

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

FCM

## Dilutions

FCM: 0.5-2 µg/Test

## Validations



Neuroligin 3 / NLGN3 Antibody, Rabbit MAb, Flow Cytometry