

GLB1L3 Polyclonal Antibody



Catalog Number: 162364

Product name

GLB1L3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to GLB1L3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 180-290 of human GLB1L3 (NP_001073876.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - GLB1L3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GLB1L3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunohistochemistry - GLB1L3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human appendix using GLB1L3 antibody at dilution of 1:150 (40x lens).