

MACC1 antibody

Catalog Number: 112364

Product name

MACC1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MACC1 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue, human colon cancer tissue.

Preparation

This antibody was obtained by immunization of MACC1 recombinant protein (Accession Number: NM_182762). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:112364(MACC1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:112364(MACC1 Antibody) at dilution of 1:50 (under 40x lens)