

Anti-TNFSF12-TNFSF13 antibody



Catalog Number: 176364

Product name

Anti-TNFSF12-TNFSF13 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to TNFSF12-TNFSF13

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tnfsf12. aa range:41-90

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of various cells using TWEAK Polyclonal Antibody diluted at 1:2000. Secondary antibody(catalog#:HA1001) was diluted at 1:20000