

# Anti-TNFSF12-TNFSF13 antibody



Catalog Number: 176364

## Product name

Rabbit IgG

Anti-TNFSF12-TNFSF13 antibody

## Applications

## Specificity

WB, IHC-P, IF

Human, Mouse

## Dilutions

## Antibody description

WB: 1:500-1:2000

Rabbit polyclonal antibody to TNFSF12-TNFSF13

IHC-P: 1:100-1:300

## Preparation

IF: 1:100-1:500

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tnfsf12. aa range:41-90

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

## Validations




Fig1:; Western Blot analysis of various cells using TWEAK Polyclonal Antibody diluted at 1:2000. Secondary antibody(catalog#:HA1001) was diluted at 1:20000