

# IL-17F antibody

Catalog Number: 111764

---

## Product name

IL-17F antibody

## Specificity

Human; other species not tested.

## Antibody description

IL-17F Rabbit Polyclonal antibody. Positive WB detected in U-937 cells, K-562 cells. Observed molecular weight by Western-blot: 32-35 kDa

## Preparation

This antibody was obtained by immunization of IL-17F recombinant protein (Accession Number: NM\_052872). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:111764(IL17F antibody) at dilution of 1:200