

# MTA3 antibody



Catalog Number: 112764

## Product name

MTA3 antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

## Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:20-1:200

## Antibody description

MTA3 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, Jurkat cells, MCF7 cells, mouse brain tissue, mouse thymus tissue, Raji cells. Positive IP detected in Raji cells. Positive IF detected in HepG2 cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 62 kDa, 68 kDa

## Preparation

This antibody was obtained by immunization of MTA3 recombinant protein (Accession Number: NM\_020744). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IF, IP


## Dilutions

Recommended Dilution:

## Validations




HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112764(MTA3 antibody) at dilution of 1:500




Immunohistochemical of paraffin-embedded human brain using Catalog No:112764(MTA3 antibody) at dilution of 1:100 (under 40x lens)

# MTA3 antibody

Catalog Number: 112764



Immunofluorescent analysis of HepG2 cells, using MTA3 antibody Catalog No:112764 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-MTA3 (IP:Catalog No:112764, 4ug; Detection:Catalog No:112764 1:500) with Raji cells lysate 2800ug.