

# ADAM5 Polyclonal Antibody



Catalog Number: 164764

## Product name

ADAM5 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ADAM5

## Preparation

Antigen: Recombinant protein of human ADAM5

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

## Validations



## Western blot - ADAM5 Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using ADAM5 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.