

# SYN2 antibody



Catalog Number: 115764

## Product name

SYN2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SYN2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 63 kDa

## Preparation

This antibody was obtained by immunization of SYN2 recombinant protein (Accession Number: NM\_133625). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115764(SYN2 Antibody) at dilution of 1:2000