

DCUN1D3 antibody



Catalog Number: 109764

Product name

DCUN1D3 antibody

Immunogen

[Human DCUN1D3 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

DCUN1D3 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, A549 cells. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of DCUN1D3 recombinant protein (Accession Number: NM_173475). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109764(DCUN1D3 Antibody) at dilution of 1:1000