

MMP7 Polyclonal Antibody



Catalog Number: 160074

Product name

MMP7 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to MMP7

Preparation

Antigen: A synthetic peptide of human MMP7

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:1000

IHC 1:50 - 1:100

IF 1:20 - 1:50

Validations



Western blot - MMP7 Polyclonal Antibody

Western blot analysis of extracts of HL-60 cells, using MMP7 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.