

C1QC Polyclonal Antibody



Catalog Number: 166074

Product name

C1QC Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to C1QC

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human C1QC (NP_001107573.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations




Western blot - C1QC Polyclonal Antibody

Western blot analysis of extracts of human serum, using C1QC antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



Immunohistochemistry - C1QC Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung using C1QC antibody at dilution of



1:100 (40x lens).


C1QC Polyclonal Antibody



Catalog Number: 166074

Immunohistochemistry - C1QC Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver using C1QC antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - C1QC Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using C1QC antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - C1QC Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using C1QC antibody at dilution of 1:100 (40x lens).