

# Anti-C Reactive Protein (1H3) Mouse antibody



Catalog Number: 168074

## Product name

Anti-C Reactive Protein (1H3) Mouse antibody

## Specificity

Human

## Antibody description

C Reactive Protein (1H3) Mouse monoclonal antibody

## Preparation

Antigen: Native C-Reactive Protein.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

## Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

ELISA

## Validations



sELISA standard Curve for C-Reactive Protein: Capture Antibody Mouse mAb to C-Reactive Protein at 1ug/ml and 168074 was used for detecting.



sELISA standard Curve for C-Reactive Protein: Capture Antibody Rabbit pAb to C-Reactive Protein at 4ug/ml and 168074 was used for detecting.