

# Human CNPY2 (His Tag) recombinant protein



Catalog Number: 500874

## General Information

### Gene Name Synonym

MIR-interacting saposin-like protein; Putative secreted protein Zsig9; Transmembrane protein 4

### Protein Construction

A DNA sequence encoding the human CNPY2 (Q9Y2B0-1) (Met1-Ser178) was expressed with a polyhistidine tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Arg 21

### Molecular Mass

The recombinant human CNPY2 consists of 173 amino acids and predicts a molecular mass of 20 KDa. It migrates as an approximately 20 KDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human CNPY2 Protein (His Tag) SDS-PAGE