

# EDN3 Polyclonal Antibody



Catalog Number: 161874

---

## Product name

EDN3 Polyclonal Antibody

## Specificity

Mouse

## Antibody description

Polyclonal antibody to EDN3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 17-238 of human EDN3 (NP\_996917.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions

IHC 1:50 - 1:200

## Validations



Immunohistochemistry - EDN3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse pancreas using EDN3 Antibody at dilution of 1:100 (40x lens).