

# GRIN2B Polyclonal Antibody



Catalog Number: 160174

## Product name

GRIN2B Polyclonal Antibody

## Specificity

Mouse

## Antibody description

Polyclonal antibody to GRIN2B

## Preparation

Antigen: A synthetic peptide of human GRIN2B

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:200 - 1:500

## Validations



Western blot - GRIN2B Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using GRIN2B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.