

LCN8 antibody

Catalog Number: 112174

Product name

LCN8 antibody

Specificity

Human; other species not tested.

Antibody description

LCN8 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of LCN8 recombinant protein (Accession Number: NM_178469). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112174(LCN8 Antibody) at dilution of 1:1000