

# Anti-RBP4 (8C12) Mouse antibody



Catalog Number: 168174

## Product name

Monoclonal

Anti-RBP4 (8C12) Mouse antibody

## Ig Type

## Specificity

Mouse IgG1

Human

## Applications

## Antibody description

ELISA

RBP4 (8C12) Mouse monoclonal antibody

## Validations

## Preparation

Antigen: Purified recombinant human RBP4 protein fragments expressed in E.coli.


## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

## Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

## Clonality



Standard Curve for RBP4: Capture Antibody Mouse mAb (8C12-C11) to RBP4 and Detector Antibody Mouse mAb(7C5-F6) to RBP4.