

ARSA antibody

Catalog Number: 108174

Product name

ARSA antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARSA Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IHC detected in human cervical cancer tissue, human testis tissue. Observed molecular weight by Western-blot: 56kd

Preparation

This antibody was obtained by immunization of ARSA recombinant protein (Accession Number: NM_001085427). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108174(ARSA antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:108174(ARSA antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:108174(ARSA antibody) at dilution of 1:50

ARSA antibody



Catalog Number: 108174

(under 40x lens)