

Anti-CD28 antibody



Catalog Number: 141574

Product name

Anti-CD28 antibody

Specificity

Human CD28/TP44

Antibody description

Rabbit Monoclonal to Human CD28

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human CD28/TP44 and conjugated with PerCP under optimum conditions, the unreacted PerCP was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of

water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



Human CD28/TP44 Flow Cytometry (FC) 15343