

# Anti-FKBP52 Rabbit antibody



Catalog Number: 167974

## Product name

Anti-FKBP52 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

FKBP52 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human FKBP52

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

WB: 1/1000

IHC: 1/200

## Validations



Western blot detection of FKBP52 in U2OS,C6 using FKBP52 Rabbit mAb(1:1000 diluted)