

# Anti-G3BP Rabbit antibody



Catalog Number: 168974

## Product name

Anti-G3BP Rabbit antibody

## Specificity

Human

## Antibody description

G3BP Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human G3BP

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/2000-1/10000

ICC/IF: 1/50-1/200


FC: 1/20

IP: 1/20

## Validations



Western blot detection of G3BP in K562, C6, HeLa cell lysates using G3BP Rabbit mAb(1:1000 diluted). Predicted band size: 52kDa. Observed band size: 68kDa.



Immunocytochemistry of G3BP(green) in HeLa cells using G3BP Rabbit mAb at dilution 1/200, and DAPI(blue)