

# CCDC41 antibody



Catalog Number: 108974

---

## Product name

CCDC41 antibody

## Immunogen

[Human CCDC41 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

CCDC41 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 82 kDa, 67 kDa

## Preparation

This antibody was obtained by immunization of CCDC41 recombinant protein (Accession Number: BC125087). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108974(CCDC41 Antibody) at dilution of 1:600