

CCDC41 antibody



Catalog Number: 108974

Product name

CCDC41 antibody

Immunogen

[Human CCDC41 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

CCDC41 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 82 kDa, 67 kDa

Preparation

This antibody was obtained by immunization of CCDC41 recombinant protein (Accession Number: BC125087). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108974(CCDC41 Antibody) at dilution of 1:600