

DOC2B antibody

Catalog Number: 109974

Product name

DOC2B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

DOC2B Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in mouse kidney tissue, human brain tissue, mouse brain tissue, human testis tissue. Observed molecular weight by Western-blot: 46kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003585).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109974(Doc2B antibody) at dilution of 1:600



Immunofluorescent analysis of HepG2 cells, using Doc2B antibody Catalog No:109974 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).