

# ELF1 antibody

Catalog Number: 110274

## Product name

ELF1 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ELF1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue, human tonsillitis tissue. Positive WB detected in PC-3 cells, A431 cells, HeLa cells, HL-60 cells, Jurkat cells, mouse thymus tissue. Observed molecular weight by Western-blot: 97 kDa

## Preparation

This antibody was obtained by immunization of ELF1 recombinant protein (Accession Number: NM\_001370332). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110274(ELF1 antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue using Catalog No:110274(ELF1 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human tonsillitis tissue using Catalog No:110274(ELF1 Antibody) at dilution of 1:50