

SHISA4 antibody



Catalog Number: 115274

Product name

SHISA4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SHISA4 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 25-27 kDa

Preparation

This antibody was obtained by immunization of SHISA4 recombinant protein (Accession Number: XM_017000419). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115274(SHISA4 antibody) at dilution of 1:2000



Immunofluorescent analysis of HEK-293 cells using Catalog No:115274(SHISA4 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG