

Anti-FOXD3 (5C4) Mouse antibody



Catalog Number: 168274

Product name

Mouse IgG1

Anti-FOXD3 (5C4) Mouse antibody

Applications

Specificity

WB

Human, Mouse, Monkey

Dilutions

Antibody description

WB: 1/500 - 1/2000

FOXD3 (5C4) Mouse monoclonal antibody

ELISA: 1/10000


Preparation

Antigen: Purified recombinant fragment of human FOXD3 expressed in E. Coli.

Validations

Formulation

Ascitic fluid containing 0.03% sodium azide.



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of FoxD3 in Jurkat and HeLa cell lysates using FoxD3 mouse mAb (1:1000 diluted). Predicted band size: 48KDa. Observed band size: 48KDa.

Clonality

Monoclonal

Ig Type